

Triple Detector Array



The Triple Detector Array employs a concentration detector, viscometer and low angle light scattering detector acting in concert, with each detector providing complimentary but different information.

The Viscotek Model 302-040 Triple Detector Array sets the standard for GPC/SEC technology and is the preferred method for characterizing natural and synthetic polymers, copolymers, proteins, excipients and other macromolecules. In a single GPC/SEC experiment, the Triple Detector Array will provide:

- Absolute molecular weight without the assumptions, extrapolations or corrections required by alternative techniques.
- Molecular size expressed as Hydrodynamic Radius (Rh) to less than 1 nm and Radius of Gyration (Rg).
- Intrinsic viscosity or molecular density.
- Information on structure, conformation, aggregation and branching.

The Model 302-040 Triple Detector Array consists of a Differential Refractive Index (RI) Detector, four-capillary Differential Viscometer and a Low Angle Light Scattering (LALS) Detector. All the detectors reside within a temperature-controlled compartment (ambient to 80°C) which also has space for 3-5 analytical GPC columns. This arrangement minimizes inter-detector volumes to reduce band-broadening effects and insures that detectors, inter-detector tubing, columns and sample reside at the same temperature throughout the course of the analysis.

for more information visit us @ www.viscotek.com

Viscotek is the leading provider of comprehensive GPC/SEC solutions, featuring multiple detection for the characterization of natural and synthetic polymers and proteins. Viscotek's innovative products and technologies are backed by expert-level technical service and analytical support.



The detectors in the Triple Detector Array are arranged in a unique series configuration as follows: light scattering – differential refractometer – differential viscometer. This configuration, with the viscometer placed after the RI, eliminates the viscometer dilution effect and break-through, hence increasing sample throughput capability and efficiency of operation.

TDA, GPCmax, FIPA, TriSEC, ARV, PolyCAL, ViscoGEL, and Viscotek are trademarks of Viscotek, Incorporated. The 4-capillary viscometer and the 2-capillary viscometer are covered by the following U.S. Patent Nos. and corresponding foreign Patent Nos.: 4,463,598; 4,627,271. Other patents pending.

©2006 Viscotek Corporation
All rights reserved.

TDA-1006

Technical Specifications

Oven

Temperature Range: Ambient–80°C

Oven Space: 7.25 inches(w) x 10.25 inches(h) x 15.25 inches(d)

Light Scattering Detector

Right Angle Light Scattering Detector uses a 90° angle geometry for maximum signal-to-noise.

Low Angle Light Scattering Detector uses a 7° angle geometry that requires no extrapolation or correction.

Cell Volume: 18 microliters

Maximum Backpressure on Cell: 15 psi

Maximum Signal: 2.5 volts

Baseline Noise: 0.4 millivolts

Baseline Drift: 3.0 millivolts/hour

Light Source: laser at 670 nm wavelength

Viscometer Detector

Patented 4 capillary, differential Wheatstone bridge configuration.

Bridge Volume: 72 microliters

Maximum Flow Rate: 1.5 ml/min (H₂O); 3.0 ml/min (THF)

DP Noise Maximum: 0.7 millivolts

DP Noise Typical: < 0.5 millivolts

DP Drift Typical: < 3.0 millivolts/hour

Sensitivity: 2.0 x 10⁻⁵ Specific Viscosity

Transducer Linearity: < 1%

Shear Rate: 3000 sec⁻¹

Refractive Index Detector

Deflection type with reference cell included.

Cell Volume: 12 microliters

Maximum Backpressure on Cell: 5 psi

Maximum Signal: 2.5 volts

Baseline Noise: 0.3 millivolts, 3.0 x 10⁻⁸ RI units

Baseline Drift: 3.0 millivolts/hour

Light Source: light emitting diode (LED) at 660 nm wavelength

Viscotek World Headquarters

15600 West Hardy Road, Houston, Texas 77060, USA

Tel: 281-445-5966 or 800-375-5966 • email: info@viscotek.com

Viscotek Europe Ltd.

20 Wellington Business Park, Crowthorne, Berkshire RG45 6LS, UK

UK +44 1344 467180 • D +49 7254 959195 • BNL +31 412 639963

F +33 4 78 51 67 08 • email: info@viscotek-europe.com