

# GPCmax™



Introducing GPCmax™, the world's first integrated pump, autosampler and degasser module specifically designed for GPC/SEC applications. It mates seamlessly with Viscotek's Triple Detector Array, Viscotek's Dual Detector\* or a simple concentration detector to form a complete GPC/SEC system.

*GPCmax™ was designed by the most experienced GPC/SEC experts in the world to deliver the performance and reliability our customers demand.*

The GPCmax™ features a variable injection volume autosampler capable of handling 2ml or 20ml vials (eliminating the need for sample transfer) and an isocratic pump for constant and pulseless solvent flow. The pump and autosampler are controlled by the Viscotek OmniSEC™ software but can be programmed manually using the front instrument panel.

The GPCmax™ also comes equipped with an in-line solvent degasser to remove dissolved gases and allow optimum performance of concentration, viscosity and light scattering detectors.

for more information visit us @ [www.viscotek.com](http://www.viscotek.com)

