

m-VROC II Chip Specification

CHIP TYPE	SAMPLE VOLUME	APPLICATIONS		VISCOSITY RANGE, mPa-s (cP)	
		VISCOSITY	SHEAR RATE	MIN	MAX
A05	Small	● ● ● ● ●	● ● ● ● ●	0.04	80
B05	Small	● ● ● ● ●	● ● ● ● ●	0.5	280
C05	Medium	● ● ● ● ●	● ● ● ● ●	3.0	1,500
C20	Large	● ● ● ● ●	● ● ● ● ●	160	100,000
E02	Medium	● ● ● ● ●	● ● ● ● ●	1.2	800
E05	Large	● ● ● ● ●	● ● ● ● ●	20	13,000